

**MALEK**  
CO<sub>2</sub> PROFESSIONAL

20  
Years of Experience

## INSPECT 3005R CO<sub>2</sub> Injector

Angiography with CO<sub>2</sub>

Picture taken with 70 ccm CO<sub>2</sub> Gas and 700 mbar pressure with 1 frame/ sec.

THE GOLD STANDARD



**MALEK**  
CO<sub>2</sub> PROFESSIONAL

Malek Medical GmbH  
Holzdamm 4  
23966 Wismar, Germany

fon +49 (0)3841-6617598  
fax +49 (0)3841-5205727  
mail@malekmedical.de

V190215

[www.malekmedical.de](http://www.malekmedical.de)



# INSPECT 3005R

## CO<sub>2</sub> Injector

Angiography with CO<sub>2</sub>

The Malek Medical Company has over 20 years experience in CO<sub>2</sub> angiography. Malek Medical is a German family company with design, development and production in Germany.

The INSPECT 3005R combines high innovation of German technology and Swiss precision.



Maximum range in pressure and volume with intelligent continuous flow over entire injection time combined with gas heating for superior results.



Pressure	Volume	Gas Heating
0 - 1400 mbar 0 - 1150 mmHg	10 - 200 ccm (ml)	36 °C

Automatic Processes: Quick Start less 45 sec  
Flushing System  
Flushing Connection Set

Cable and Controls: Power Cable up to 8m  
Foot Switch Cable up to 8m  
Remote Control optional

### INSPECT 3005R - The Gold Standard

Superior results in upper extremities and infradiaphragmal. By the inventor of DSA with CO<sub>2</sub>.

CO<sub>2</sub> Injectors developed by CEO Frank Malek are worldwide in over 100 hospitals state of the art for more than 20 years.

### Head of Vascular Surgery, Germany

„CO<sub>2</sub> Angiography with the CO<sub>2</sub> Injector by Malek Medical is established in our clinic as «First-choice»-standard for all intraoperative, endovascular procedures since already many years. A comfortable handling, interconnection with X-ray unit and preheating of the gas ensure a obviously increased safety for patients and doctors.“

